
INQUA® Vernebler

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REF BR021200

CE 0123

| | | |
|-----------|-----------------------------|------------|
| de | Gebrauchsanweisung | 3 |
| en | Instructions for use | 29 |
| es | Instrucciones de uso | 55 |
| fr | Notice d'utilisation | 81 |
| it | Istruzioni per l'uso | 107 |
| nl | Gebruiksaanwijzing | 133 |
| pt | Manual de Instruções | 159 |
| tr | Kullanım talimatları | 185 |

INQUA® Vernebler

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Instructions for use

TABLE OF CONTENTS

| | |
|-------------------------------------------------------------------|-----------|
| GENERAL | 30 |
| Manufacturer's data | 30 |
| Notes on the instructions for use | 30 |
| PRODUCT DESCRIPTION | 32 |
| Components | 32 |
| Correct use | 33 |
| Conformity | 34 |
| SAFETY INSTRUCTIONS | 35 |
| APPLICATION | 36 |
| Preparing for inhalation | 37 |
| Performing the inhalation | 39 |
| Completing the inhalation | 42 |
| HYGIENIC RE-USE | 43 |
| Care of the connection tubing | 44 |
| Preparation for cleaning/disinfection | 45 |
| Cleaning with warm drinking water | 45 |
| Disinfection | 46 |
| Drying | 48 |
| Storage | 48 |
| MISCELLANEOUS | 49 |
| Storage and transportation conditions | 49 |
| Material information | 49 |
| Operating life | 49 |
| Disposal | 50 |
| Legend | 50 |
| SPARE PARTS/ACCESSORIES | 51 |
| TECHNICAL DATA | 51 |
| General | 51 |
| Aerosol characteristics in accordance with EN 13544-1 appendix CC | 51 |

GENERAL

WARNING:

- Read these instructions carefully and completely before using for the first time, and keep them in a safe place. Do not use the INQUA® nebuliser unless you fully understand the contents of these instructions for use.
- If symptoms persist or if your condition worsens, discontinue the application and contact your doctor immediately.

Manufacturer's data

See back of instructions for use.

Notes on the instructions for use

Purpose of the instructions for use

These instructions for use describe how the INQUA® nebuliser is assembled, how it works, and how it should be used. They are intended to ensure that the individual components can be used in complete safety.

Instructions for using the compressor that will be attached to the nebuliser are not included in these instructions for use. Please read the instructions for use of the compressor which accompany that item. If necessary, you can order another copy of these instructions for use from INQUA GmbH (see manufacturer information on the back of the instructions for use).

Representation of hazard warnings

In these instructions for use, hazard warnings are indicated by a warning triangle and a signal word. The following signal words are used depending on the nature and degree of danger being described.

- **NOTICE**

Indicates general precautionary measures which are to be observed to avoid damaging the product during use.

- **CAUTION**

Indicates a hazard that may cause mild to moderate injury or damage if it is not avoided.

- **WARNING**

Indicates a hazard that may cause serious injury or death if it is not avoided.

Forms of emphasis


The following forms of emphasis are used in these instructions for use to indicate important parts of the text:

- ⇒ The arrow indicates that a given action must be taken.

PRODUCT DESCRIPTION

Components

INQUA® nebuliser set (applied part)

| | | | |
|----------|------------------------|----------|---------------------------------------------------------|
| 1 | Nebuliser (pre-fitted) | 5 | Connection tubing |
| 2 | Mouthpiece | 6 | Air filter for the compressor manufactured by INQUA® |
| 3 | Adult mask | 7 | Removal tool for changing the filter |
| 4 | Child mask | | |

Correct use

Together with the compressor being used, the INQUA® nebuliser forms an inhalation device for treatment of the upper and lower respiratory tracts with nebulised medication. It is suitable for temporary¹⁾ oral or nasal inhalation. It may be used to inhale all medications that have been authorised for use in inhalation therapy or prescribed or recommended by a doctor. Always read the notes (restrictions, frequency of application) in the accompanying information for the medication to be inhaled, and ask your doctor or pharmacist if you have any questions.

The INQUA® nebuliser generates a moist-cool aerosol that cannot be heated up. It must be used exclusively for inhalation treatment, not in conjunction with anaesthesia breathing systems or in the breathing system of respirators. Any application other than the use described here is considered improper, and therefore dangerous. In all cases the user is liable for the safe operation of the nebuliser if it is not used according to the instructions.

The INQUA® nebuliser is intended exclusively for use in the home and together with the INQUA® compressor. For reasons of hygiene, it must only be used by one person.

Note that the connection tubing included with the nebuliser kit must only be used for the inhalation device. If several tubing systems are used for different applications at the same time (e.g. infusion and inhalation), there is a danger that the tubing systems may be confused. In extreme cases, this can even result in death.

The nebuliser is only designed for patients who are able to breathe by themselves and are conscious. This is the only way to ensure effective treatment and avoid a choking hazard.

1) Each application lasts less than 1 hour according to the European Medical Device Directive 93/42/EEC, in this case max. 20 minutes.

The nebuliser system is **not** suitable for patients who have undergone a tracheostomy. Patients who have undergone a tracheostomy breathe through a surgical opening in the windpipe (tracheostoma). It is not possible to assure effective treatment for this clinical picture with this equipment.

The nebuliser equipment is suitable for use by children aged three and older up to adult. The nebuliser system can be used for inhalation treatment of the **lower respiratory tract** (bronchi) if the patient uses the **mouthpiece**, breathing in through the mouth and out through the nose. The mask should only be used for treatment of the lower respiratory tract in exceptional circumstances. In this case, the patient must be sure to breathe in through the mouth.

The nebuliser system can be used for inhalation treatment of the **upper respiratory tract** if the patient is a child using the accompanying **child mask** and breathing in and out through the nose. Adults use the **nasal joining piece**, breathing in and out through the nose.

Inhalation treatment for babies: In order to ensure efficient inhalation treatment for the lower respiratory tract in babies (up to age three), you need a special baby nebuliser. In this case, contact the manufacturer (for address see the back of the instructions for use).

Conformity

The INQUA® nebuliser satisfies the requirements of Directive 93/42/EEC (EC Medical Device Directive).

SAFETY INSTRUCTIONS

The following general instructions apply:

- People with mental, physical or sensory disabilities and children do not recognise the risks that may arise when using the inhalation device and its components. Never allow anyone from these groups access to the inhalation device without supervision.

For example, there is a risk that individual parts may be swallowed and block the respiratory tract, or bitten or bitten off, or that the inhalation device may be pulled off the table by the connection tubing or power cord. There is also a danger of strangulation with the connection tubing or power cord. Therefore, make sure that you always keep the compressor, the nebuliser and the accessories out of the reach of babies and infants.

- INQUA GmbH cannot be held responsible for damage that is caused by improper or incorrect use, e.g., use of a compressor made by another manufacturer. Therefore, you are advised to read the instructions for use of the INQUA® nebuliser and the INQUA® compressor carefully and to familiarise yourself with the device before using it. **It is particularly important to follow the instructions for use of the INQUA® compressor.**
- You must conscientiously follow the general hygiene procedures when using your nebuliser, e.g., thorough washing of hands or disinfection if required.
- Use only liquid medications that have been approved for use in inhalation therapy with a compressed air inhalation device in your INQUA® nebuliser. Ask your doctor or pharmacist about the selection and application of the medications available.
- In general, the following are not suitable for use in inhalation devices: cough mixtures, gargles, balsamic preparations, oils of medicinal plants, and drops that are used as ointments or added to the bath. Inhaling these medications can cause irreversible damage to the respiratory tracts.
Such substances can also impair the function of the inhalation device permanently.

- Use only accessories and spare parts that are recommended by INQUA in order to obtain the best possible treatment (see Chapter "Spare parts/Accessories", page 51).
- It is most important to follow the instructions in chapter "Hygienic re-use", page 43.
- Note the required ambient conditions for operation:

| | |
|-----------------------|-----------------------------|
| Ambient temperature: | +10 °C to +40 °C |
| Relative humidity: | 30% to 75% (non-condensing) |
| Atmospheric pressure: | 700 hPa to 1,060 hPa |

APPLICATION

WARNING: Danger of infection from contaminated nebuliser

Follow the hygiene regulations before each use and after any prolonged interruption in use. Make sure that the nebuliser has been cleaned and disinfected before it is used for the first time or after any previous treatment (see Chapter "Hygienic re-use", page 43).


CAUTION: Impaired treatment due to damaged parts

Check all nebuliser components before each use. Replace any broken, misshapen or seriously discoloured parts. Please also follow the instructions for assembly given below. Damaged components and/or an incorrectly assembled nebuliser may impair **functioning of the nebuliser** and thus **treatment as well**.

Preparing for inhalation

Assembling the nebuliser

If the nebuliser has been dismantled, follow the instructions below to assemble it for use:


Note that the nebuliser upper section (**1a**) and the nebuliser lower section (**1b**) are attached to each other by a "bayonet" closure.

- ⇒ The nebuliser upper section is inserted carefully into the nebuliser lower section.
- ⇒ Pay particular attention to the indication **front** in the drawing.

The lugs in the nebuliser upper section are different sizes and only fit in the different sized slots in the nebuliser lower section in a certain direction.

Twist the nebuliser upper section clockwise while pressing gently to engage the lugs on the nebuliser upper section in the slots in the lower section. The two nebuliser sections will be attached securely by gently twisting clockwise.

⇒ Connect one end of the connection tubing to the air outlet on the INQUA® compressor (see the instructions for use of the INQUA® compressor) and the other end to the INQUA® nebuliser.

Inhalation for the lower respiratory tract:

⇒ Attach the mouthpiece or child mask to the nebuliser. The child mask should only be used in exceptional cases.

Inhalation for the upper respiratory tract:

⇒ Attach the mask to the nebuliser.

Filling with inhalation solution

⇒ Insert the nebuliser in the holder on the compressor intended for this purpose.
⇒ Pour the quantity of inhalation solution prescribed by your doctor or the manufacturer into the top of the nebuliser. Make sure that the liquid added does not pass the top scale graduation mark (max. fill quantity 8 ml). Overfilling can greatly impair the nebulisation of medication.

The nebuliser is ready for use.

Performing the inhalation

WARNING: Danger to health from improper use

Small parts can be swallowed and block the respiratory tract, thus posing a choking hazard.

Please note that infants and anyone who requires assistance must be supervised constantly during the inhalation treatment.

CAUTION: Inadequate medication nebulising due to incorrect assembly

Before the patient starts inhaling, ensure that all parts are firmly connected to each other. Make sure that the connection tubing is fitted securely on the air outlet on the INQUA® compressor and on the nebuliser.

CAUTION: Inadequate delivery of medication because nebuliser held incorrectly

Always hold the nebuliser upright during inhalation.

Inhalation for the lower respiratory tract

Inhalation with the mouthpiece is the most effective form of inhalation for the lower respiratory tract, since it results in minimum medication loss on the way to the lung. The child mask is only recommended for children if they are unable to inhale using the mouthpiece. When a mask is used, the patient must always breathe in through the mouth during the session.

- ⇒ Sit in an upright position and relax.
- ⇒ Switch the compressor on.
The nebuliser will generate a visible mist (aerosol).

- ⇒ Hold the mouthpiece between your teeth and enclose it with your lips.

- ⇒ Breathe in through the mouthpiece as calmly, deeply and evenly as possible.

- ⇒ After a few seconds, breathe out through your nose.
When breathing out through the nose, a small quantity of mist (aerosol) may escape from the top of the nebuliser. These losses are not significant.

Inhalation for the upper respiratory tract

- ⇒ Hold the nebuliser vertically, sit in a relaxed, upright position.
- ⇒ Switch the compressor on.
The nebuliser will generate a visible mist (aerosol).

When using the mask:

- ⇒ Hold the mask up to the face so that it forms an airtight seal over the mouth and nose.
- ⇒ Breathe in calmly, deeply and evenly through the nose, and after a few seconds breathe out again through the nose.

Completing the inhalation

The duration of inhalation depends on the quantity of medication used and the breathing pattern.

Inhalation is complete when no more aerosol is produced by the nebuliser, or the aerosol only emerges from the nebuliser intermittently. The noise made by the nebuliser will change.

A small quantity of medication (approx. 1 ml) will be left in the nebuliser even after an inhalation session.

- ⇒ When the inhalation session has ended, turn the INQUA® compressor off at the On/Off switch.
- ⇒ Pull the power plug out of the socket (see instructions for use of the INQUA® compressor) and disconnect the nebuliser from the connection tubing.

If the nebuliser is to be used for several different medications one after the other in one inhalation session, the nebuliser must be cleaned before it is filled with different medication (see Chapter "Hygienic re-use", Section "Cleaning with warm drinking water", page 45).

WARNING: Bacterial contamination caused by residues

Bacteria can breed very quickly in nebulisers, mouthpieces or masks that have not been cleaned properly and in liquids that are left in the nebuliser, thereby endangering health, for example by causing infections.

You should therefore begin preparation for hygienic re-use as soon as the inhalation session is ended. Always empty out the residue in the nebuliser and never store any liquids in it.

HYGIENIC RE-USE

WARNING:

In order to prevent a risk to health, e.g., infection due to a contaminated nebuliser, it is imperative that the following hygiene regulations are followed.

The INQUA® nebuliser is designed for multiple uses. It must not be used by more than one person. A separate nebuliser must be provided for each patient.

To do this, follow these instructions:

- The nebuliser, the mouthpiece and the child or adult mask must be cleaned **and** disinfected before they are used for the first time and after each application.
- Wash your hands thoroughly before every use and cleaning. Additional requirements concerning the necessary hygienic preparations (hand disinfection, handling of medications and inhalation solutions) in high-risk groups (e.g., users with suppressed immune systems, patients with cystic fibrosis) can be obtained from the respective self-help groups.
- The connection tubing cannot be cleaned nor disinfected either chemically or by boiling (see Chapter "Care of the connection tubing", page 44).
- Ensure that the INQUA® nebuliser is dried thoroughly after each cleaning and disinfection. Condensation or residual wetness can present an increased risk through the growth of bacteria.
- Please read Section "Material information" in Chapter "Miscellaneous" on page 49.
- Check the parts of your nebuliser regularly and replace any defective (broken, misshapen, discoloured) parts (see Chapter "Spare parts/Accessories", page 51).
- Please also read and follow the hygiene instructions in the instructions for use of the INQUA® compressor.

Like all plastic parts, the INQUA® nebuliser is susceptible to a certain degree of wear when it is used and reprocessed frequently (wearing part). Over time, this can cause a change in the aerosol, which in turn impairs the therapeutic efficiency.

You should replace the nebuliser at least once a year. Replacement sets can be obtained from your sales partner or from INQUA GmbH (see Chapter "Spare parts/Accessories", page 51).

Proof of suitability of the INQUA® nebuliser for effective cleaning and disinfection has been obtained through an independent testing laboratory using the procedures recommended starting on page 45.

Using the named alternatives remains the responsibility of the user.

Care of the connection tubing

- ⇒ After inhalation, check the connection tubing for condensation.
To do this, after each treatment pull the tubing off of the nebuliser and allow the INQUA® compressor to continue running until all moisture has been removed by the airstream.
- ⇒ Do not use water or any other cleaning or disinfecting solutions to clean the connection tubing.
Replace the tubing if it is contaminated.


WARNING: Bacterial contamination caused by residual moisture

Any residual moisture in the connection tubing can present a risk to health, e.g., by causing an infection.

Preparation for cleaning/disinfection

Directly after each treatment, all parts of the INQUA® nebuliser must be cleaned of medication residue and contamination.

- ⇒ If this has not been done, disconnect the connection tubing from the INQUA® nebuliser.
- ⇒ Dismantle the nebuliser into its individual parts.

In order to help reduce water pollution, do not dispose of the medication residue in the waste water system. Wipe it up with tissues and dispose of the tissues with non-recyclable waste.

Cleaning with warm drinking water

⇒ Clean all parts of the nebuliser thoroughly for 5 minutes with warm drinking water (approx. 40 °C) and a little dishwashing liquid with no skin protecting or disinfecting additives (dosage according to the manufacturer's recommendations) until all residues of the inhalation solution and other contaminants have been removed.

If necessary, use a small brush to remove dirt. For hygiene reasons, this brush must then be used exclusively for cleaning the nebuliser.

- ⇒ Then rinse all parts thoroughly (for at least 1 minute) under warm running water (approx. 40 °C without dishwashing liquid) to remove all remaining dishwashing liquid.
- ⇒ Pour out the water that has collected in the individual parts. You can remove excess water more quickly by shaking the parts.

Disinfection

Disinfect the dismantled nebuliser and the mouthpiece or the mask used immediately after cleaning.

Disinfection with boiling water

- ⇒ Place the individual parts of your INQUA® nebuliser in boiling water for at least 15 minutes.

Use a clean pot and fresh drinking water.

Drinking water contains dissolved particles (e.g., salts) that can form a milky film on plastic parts during the boiling process, depending on their type and quantity. This film can be prevented by using distilled water in the boiling operation.

When boiling, make sure that there is enough water in the pot to prevent the nebuliser parts from coming into contact with the base of the pot.

WARNING: Risk of fire if boiling operation is not supervised!

As soon as the water has evaporated completely, the plastic in the pot will start to melt. If the heat source is not cut off, the plastic will decompose, giving off smoke, and may cause a fire.

Never leave the boiling procedure unattended!

Possible alternative: Disinfection with hot steam

The INQUA® nebuliser can be disinfected using a steam disinfecter (disinfecter for baby milk bottles). To ensure that the hot steam passes through all the individual components of the nebuliser, the components must be placed upright in the disinfecter. For effective disinfection, use a disinfecter with a runtime of at least 6 minutes. Regarding disinfection, the duration of the disinfection procedure and the quantity of water required for this, follow the instructions for use of the disinfecter you are using.

Microwave disinfectors are not suitable for disinfecting the nebuliser and its components (mouthpiece, masks).

NOTICE:

Microwave disinfectors are not suitable for disinfecting the nebuliser and its components (mouthpiece, masks).

CAUTION: Risk of infection due to inadequate disinfection

Inadequate disinfection encourages the growth of bacteria and thus increases the risk of infection. Thorough disinfection has not been completed until the disinfecter automatically switches itself off or the minimum disinfection time specified in the instructions for use of the disinfecter has elapsed. Therefore, do not switch the device off prematurely. Also make sure that the disinfecter is kept clean and check that it is working perfectly regularly.

Drying

- ⇒ After disinfection, place the INQUA® nebuliser components on a dry, clean and absorbent surface.
You can remove excess water more quickly by shaking all parts.
- ⇒ Let them dry completely.

The components should not be dried in damp rooms (e.g., the bathroom).

Storage

WARNING: Bacterial contamination caused by residual moisture

Inadequate drying can present a risk to health, e.g. by causing an infection.

Make sure that all components are free of residual moisture before reassembling them and/or packing them away or using them for inhalation again. All components must be perfectly dry.

- ⇒ Check all product components. Replace any broken, misshapen or discoloured parts.
- ⇒ Between uses, and in particular during long breaks in treatment, wrap the dry nebuliser in a dry, clean, lint-free cloth (e.g. a tea-towel) and keep it in a dry, dust-free environment.

WARNING: Danger to health from choking or misuse

If children operate the inhalation device without supervision, there is a risk that they may swallow small parts, which may block the respiratory tract, or that incorrect use may cause injury (e.g. strangulation with the connection tubing or power cord).

Keep the INQUA® nebuliser out of children's reach.

MISCELLANEOUS

Storage and transportation conditions

Do not store the INQUA® nebuliser in damp rooms (e.g. bathroom).

Do not store or transport it together with damp objects.

Protect the INQUA® nebuliser from continuous, direct sunlight while it is being stored and moved.

In general, the nebuliser should be stored as described in Chapter "Hygienic re-use", Section "Storage", page 48 between applications. Be sure to observe the following ambient conditions:

| | |
|------------------------------|--------------------------------------------------------------|
| Minimum ambient temperature: | -25 °C (without monitoring of relative humidity) |
| Maximum ambient temperature: | +70 °C (with relative humidity of up to 93%, non-condensing) |
| Humidity: | 0% – 93% |
| Atmospheric pressure: | 500 hPa to 1,060 hPa |

Material information

The nebuliser upper part, the nebuliser lower part, the masks and the mouthpiece are made from polypropylene.

Operating life

The INQUA® nebuliser has an average expected life cycle of:

- approx. 300 disinfection cycles
- max. 1 year

Disposal

The components of the INQUA® nebuliser and the replacement set must be disposed of with non-recyclable waste.

Legend

Symbol for manufacturer

Symbol for Follow instructions for use

The INQUA® nebuliser satisfies the basic requirements as set forth in Appendix I of Directive 93/42/EEC concerning medical devices.

Symbol for reference number

Symbol for batch number

Symbol for Caution

SPARE PARTS/ACCESSORIES

Applied part

| Item description | Item no. |
|-------------------------------------------|----------|
| INQUA® nebuliser set for INQUA® Inhalator | BR021200 |

TECHNICAL DATA

General

| | |
|-----------------------------------|-------------------|
| Recommended nebuliser gases: | Air |
| Minimum flow: | (3.0 l/min) |
| Corresponding operating pressure: | 50 kPa (0.5 bar) |
| Maximum flow: | (6.0 l/min) |
| Corresponding operating pressure: | 200 kPa (2.0 bar) |
| Minimum fill volume: | 2 ml |
| Maximum fill volume: | 8 ml |


Aerosol characteristics in accordance with EN 13544-1 appendix CC

The output of the nebuliser depends on the medication being used. Consequently, the following data may not apply for medications that exist as suspensions or are very viscous. If in doubt, ask the manufacturer of the medication concerned for advice.

Determination of particle size

Fill the INQUA® nebuliser with 3 ml of a 2.5% test solution and then attach the INQUA® nebuliser to an impactor (NGI™ = Next Generation Impactor™) and run for 3 minutes.

Cumulative mass distribution for the INQUA® nebuliser at the minimum flow (3.0 l/min) and the maximum flow (6.0 l/min):

(Averages from each of 2 separate measurement series with 3 nebuliser/compressed air combinations)

Determination of aerosol output¹⁾

Fill the nebuliser with 2 ml of a 1% test solution. Attach the nebuliser to the Breath Simulator until the aerosol becomes irregular, then one more minute.

Aerosol output in [ml] for minimum flow (3.0 l/min) and maximum flow (6.0 l/min):

| | |
|-----------------------|---------|
| Min. flow (3.0 l/min) | 0.21 ml |
| Max. flow (6.0 l/min) | 0.19 ml |

(Averages from each of 2 separate measurement series with 3 nebuliser/compressed air combinations)

1) Aerosol quantity output by the nebuliser system for a specified fill volume

Determination of the aerosol output rate¹⁾

Fill the nebuliser with 2 ml of a 1% test solution. Attach the nebuliser to the Breath Simulator for at least 1 minute.

Aerosol output rate in [ml/min] for minimum flow (3.0 l/min) and maximum flow (6.0 l/min):

| | |
|-----------------------|--------------|
| Min. flow (3.0 l/min) | 0.096 ml/min |
| Max. flow (6.0 l/min) | 0.144 ml/min |

(Averages from each of 2 separate measurement series with 3 nebuliser/compressed air combinations)

We reserve the right to make technical changes.

Information as of: October 2013

1) Quantity of aerosol that is output by the nebuliser system per unit of time

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