
# INQUA® Inhalator

REF BR021000

CE 0123

<b>de</b>	<b>GEBRAUCHSANWEISUNG</b>	<b>5</b>
<b>en</b>	<b>INSTRUCTIONS FOR USE</b>	<b>27</b>
<b>es</b>	<b>INSTRUCCIONES DE USO</b>	<b>49</b>
<b>fr</b>	<b>NOTICE D'UTILISATION</b>	<b>75</b>
<b>it</b>	<b>ISTRUZIONI PER L'USO</b>	<b>99</b>
<b>nl</b>	<b>GEBRUIKSAANWIJZING</b>	<b>121</b>
<b>pt</b>	<b>MANUAL DE INSTRUÇÕES</b>	<b>145</b>
<b>tr</b>	<b>KULLANIM TALİMATLARI</b>	<b>169</b>

# **INQUA® Inhalator**

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## **Instructions for use**

### **TABLE OF CONTENTS**

<b>1</b>	<b>GENERAL</b>	<b>28</b>
1.1	Manufacturer's data .....	28
1.2	Notes on the instructions for use .....	28
<b>2</b>	<b>PRODUCT DESCRIPTION</b>	<b>30</b>
2.1	Notes on the product .....	30
2.2	Technical data .....	32
2.3	Components .....	34
<b>3</b>	<b>NOTES ON PRODUCT SAFETY</b>	<b>35</b>
<b>4</b>	<b>APPLICATION</b>	<b>37</b>
4.1	Preparation .....	37
4.2	Installation location .....	37
4.3	Preparing for inhalation .....	40
4.4	Performing the inhalation.....	41
4.5	Completing the inhalation .....	42
<b>5</b>	<b>CLEANING, MAINTENANCE, REPAIR</b>	<b>43</b>
5.1	Cleaning tasks .....	43
5.2	Maintenance .....	44
5.3	Repairs .....	45
5.4	Storage .....	45
<b>6</b>	<b>TROUBLESHOOTING</b>	<b>45</b>
<b>7</b>	<b>MISCELLANEOUS</b>	<b>46</b>
7.1	Transportation .....	46
7.2	Storage .....	46
7.3	Terms and conditions of warranty .....	46
7.4	Disposal .....	48
7.5	Information about RoHS .....	48
<b>8</b>	<b>SPARE PARTS</b>	<b>48</b>

# 1 GENERAL

Read these instructions for use carefully and completely before using for the first time. Do not use the INQUA® Inhaler until you have read and fully understood the contents of these instructions for use and the contents of the instructions for use of the INQUA® nebuliser. If symptoms persist or if your condition worsens, discontinue the application and contact your doctor immediately.

---

## 1.1 Manufacturer's data

See back of instructions for use.

## 1.2 Notes on the instructions for use

### Purpose of the instructions for use

These instructions for use describe how the INQUA® Inhalator is assembled, how it works, and how it should be used. They are intended to ensure that the individual components can be used in complete safety.

These instructions for use do not include information on using the nebuliser that is to be attached to the compressor, nor for the masks. Please read the instructions for use of the nebuliser which accompany that item.

### Forms of emphasis

The following forms of emphasis are used in these instructions for use to indicate important parts of the text:

- ⇒ The arrow indicates that a given action must be taken.

### Representation of hazard warnings

In these instructions for use, hazard warnings are indicated by a warning triangle and a signal word. The following signal words are used depending on the nature and degree of danger being described.

- **WARNING**

Indicates a hazard that may cause serious injury or death if it is not avoided.

- **CAUTION**

Indicates a hazard that may cause mild to moderate injury or damage if it is not avoided.

- **NOTICE**

Indicates general precautionary measures which are to be observed to avoid damaging the product during use.


### **Note on storing the instructions for use**

These instructions for use must be stored so that anyone who is using the INQUA® Inhalator can access them whenever necessary. They are an integral part of the product and must be included with the device if the INQUA® Inhalator changes hands (e.g. is sold). The instructions for use are available in printed form on request from INQUA GmbH. **Check whether a more recent version of the instructions for use is available at regular intervals.**

## 2 PRODUCT DESCRIPTION

### 2.1 Notes on the product      Product components


Figure 1

- |   |                          |   |                   |
|---|--------------------------|---|-------------------|
| 1 | INQUA® compressor        | 3 | Connection tubing |
| 2 | Nebuliser (applied part) |   |                   |

### Control elements

Figure 2

- |   |   |    |                      |
|---|---|----|----------------------|
| 4 | Power cord                                      | 8  | On/off switch        |
| 5 | Carrying handle                                 | 9  | Holder for nebuliser |
| 6 | Air filter                                      | 10 | Storage compartment  |
| 7 | Air connection (marked "Air" on the compressor) |    |                      |

## **Intended use**

The purpose of the compressor is to generate compressed air for operating an INQUA nebuliser.

## **Correct use**

When used in conjunction with the accompanying nebuliser, the INQUA® compressor forms an inhalation device for treatment of the upper and lower respiratory tracts with nebulised medication. It is designed for temporary<sup>1)</sup> oral or nasal inhalation. It can be used to inhale all liquid medications which have been approved for use in aerosol therapy or which have been prescribed or recommended by a qualified physician. Always read the notes (restrictions, frequency of application) in the accompanying information for the medication to be inhaled, and ask your doctor or pharmacist if you have any questions.

The INQUA® Inhalator is intended exclusively for use in the home. It can be operated by the patient himself.

The INQUA® Inhalator must only be used for inhalation treatment. Any application other than the use described here is considered improper, and therefore dangerous.

In all cases the user is liable for the safe operation of the appliance if it is not used according to the instructions.

## **Contraindications**

None known.

## **Conformity**

The INQUA® Inhalator satisfies the requirements of Directive 93/42/EEC (EC Medical Device Directive).

## **Marking on the product**


The following markings are on the compressor:


Please follow the instructions for use  
(background = blue, symbol = white)

Type identification plate  
(example, bottom of compressor)

"Off" for a part of the device

"On"

1. Each application lasts less than 1 hour according to the European Medical Device Directive 93/42/EEC, in this case max. 20 minutes.

---

~	Alternating current
	Protection class of the part used: Type BF
	Protect the device from moisture
	Protection class II device
	Minimum and maximum ambient temperature
	Minimum and maximum humidity
	Minimum and maximum air pressure
<b>REF</b>	Order number
<b>SN</b>	Serial number
	The medical device was distributed commercially after 13 August 2005. <b>This product must not be disposed of with normal domestic waste.</b> The symbol of the refuse bin with a cross through it indicates that it must be collected separately.
	Manufacturer
	This product satisfies the requirements of 93/42/EEC (Medical devices) and 2011/65/EU (RoHS).

---

## 2.2 Technical data General

Dimensions:	L 28.5 cm, W 22.5 cm, H 13.0 cm
Weight:	2.10 kg
Max. sound intensity:	59 dBA acc. to: DIN EN 13544-1, Section 26

### Electrical data

Supply voltage:	230 V
Mains frequency:	50 Hz
Power consumption:	0.7 A

### Product specifications

Jet flow:	4 l/min
Operating pressure:	1.1 bar

## **Electromagnetic compatibility**

The INQUA® Inhalator satisfies the requirements of DIN EN 60601-1-2 with regard to electromagnetic compatibility.

## **Safety mechanisms**

The compressor is equipped with a fuse and a temperature monitor to prevent the compressor from overheating.

## **Operating life**

The compressor has an average expected operating life of:

- approx. 1000 operating hours
- up to 4 years

## **Environmental conditions for operation**

---

Ambient temperature:	+10 °C to +30 °C
----------------------	------------------

---

Relative humidity:	30% to 75% (non-condensing)
--------------------	-----------------------------

---

Atmospheric pressure:	700 hPa to 1060 hPa
-----------------------	---------------------

---

Use of the device in areas with elevated magnetic or electrical radiation (e.g., close to an MRI scanner) is not permitted. The device is intended for operation in healthcare situations at home. It is not permitted to operate it in trains, motor vehicles or aeroplanes.

## **Classification of the INQUA® Inhaler according to DIN EN 60601-1**

---

Type of electric shock protection:	Protection class II
Degree of protection from electric shock of the applied part:	Type BF
Degree of protection against water ingress in accordance with IEC 529 (IP rating):	No protection
Degree of protection when used in the presence of flammable mixtures of anaesthetics with air or with oxygen or nitrous oxide:	No protection
Operating mode:	Continuous operation

---

We reserve the right to make technical changes.

## **2.3 Components**

The INQUA® Inhalator item no. BR021000 will be delivered with the following parts:

<b>Number</b>	<b>Description</b>
1	INQUA® compressor
1	Nebuliser kit (applied part), consisting of:
	1 Nebuliser
	1 Mouthpiece
	1 Adult mask
	1 Child mask
	1 Connection tubing
	1 Air filter for the compressor
	1 Removal tool for changing the filter

### **Upon receipt**

- ⇒ Check the components.
- ⇒ Inspect the individual parts: Has the compressor suffered visible damage during shipping? If the compressor is visibly damaged or if you suspect that there is a fault, do not operate the device under any circumstances. Notify the manufacturer (see back of instructions for use).
- ⇒ Keep the delivery note and the original packaging in a safe place.

### **3 NOTES ON PRODUCT SAFETY**

During development of the INQUA® Inhaler, great importance was attached to compliance with the recognised safety regulations and reflecting the latest state of the art. However, improper use of the product may endanger the user or the patient. In order to avoid this, follow the instructions in the following chapters.

#### **The following general instructions apply:**

- People with mental, physical or sensory disabilities and children do not recognise the risks that may arise when using the inhalation device and its components. Never allow anyone from these groups access to the inhalation device without supervision. There is a risk for example that individual parts may be swallowed, chewed or bitten off, or that the inhalation device is pulled off a table by the connection tube or the power cord. There is also a risk of strangulation with the connection tube or the power cord.
- Small parts such as the air filter or the removal tool can be swallowed or block the airways and so pose a danger of suffocation. Therefore, keep and dispose of all small parts out of the reach of children!
- Only plug the compressor into mains sockets that are compliant with current safety regulations. If necessary, consult a qualified electrician. If the compressor is plugged into a non-compliant mains socket, there is a danger of electric shock.
- The compressor is designed for use with the supply voltage printed on the type identification plate. If the compressor is plugged into a mains socket with a different supply voltage, the compressor may be damaged.
- To avoid the risk of electric shock, disconnect the power plug from the mains socket immediately
  - if a malfunction occurs during operation,
  - before cleaning and maintaining the appliance,
  - after each use.
- Do not pull the power cord over sharp edges. Do not allow it to become pinched or hang unsupported. There is a risk of electric shock if the compressor is operated with a damaged power cord.
- Do not attempt to make any conversions or unauthorised additions to the product, because this will impair the product's operating reliability.
- Store and operate the inhalation device only in places where it is supported firmly and cannot fall. Because of its weight, a falling inhalation device poses a danger for any person, animal or object underneath it.

- Use only accessories recommended by INQUA to operate the INQUA® Inhalator. Operational reliability and the success of the treatment may be compromised if accessories from other manufacturers are used.
- Never leave the compressor unsupervised while it is in operation.
- The compressor is not protected against sprayed water. Protect the compressor, the power cord and the power plug from liquids and moisture.
- **Please follow the instructions for use of the accompanying nebuliser as well!**

## 4 APPLICATION

### 4.1 Preparation

#### **WARNING: Danger of electric shock if compressor is damaged**

The compressor has been designed so that the user is protected from live parts. This protection may not be assured in the case of a defective compressor.

- ⇒ If the compressor, the power cord or the power plug is visibly damaged or if you suspect a fault, do not operate the compressor under any circumstances. Disconnect the power plug from the mains socket immediately and notify the manufacturer immediately (see back of instructions for use).
- ⇒ Check that the compressor housing, the power cord and the power plug are in perfect condition **before** each inhalation session.
- ⇒ Check the connection tubing and the air filter for dirt, discolouration or dampness. Faulty, damaged or dirty parts must be replaced before the device is used.

### 4.2 Installation location

#### **CAUTION: Fire hazard**

The compressor is an electrical device. Damaged electrical components can cause a fire. To avoid the risk of fire, follow the instructions below:

- ⇒ Do not operate the compressor in explosion-prone areas or in the presence of oxidising gases (e.g. oxygen, nitrous oxide, flammable anaesthetics), or in damp locations.
- ⇒ Protect the power cord from excessive heat and oil.
- ⇒ Do not operate the compressor close to flammable objects such as curtains, tablecloths, or paper.
- ⇒ Keep the compressor and the power cord away from hot surfaces (e.g., cookers, radiators, open fires). These might damage the compressor housing or the power cord insulation.
- ⇒ Keep the power cord away from household pets (e.g. rodents). Animals may damage the insulation on the power cord.
- ⇒ Plug the compressor into an easily accessible mains socket so that it can be disconnected from the power source quickly.

#### **CAUTION:Risk of injury from tripping or caused by the device falling**

A compressor placed in a poor position or for which the cables are routed untidily poses a risk of injury.

- ⇒ Route the power cord in such a way that it does not present a tripping hazard and people cannot become entangled in it.
- ⇒ Make sure that children are not able to pull the compressor down by the power cord or connection tubing.

#### **CAUTION:Quality of treatment impaired by electromagnetic interference**

Electrical devices can cause electromagnetic interference. Interference can impair the function of the devices and thus also the effectiveness of the treatment.

- ⇒ Do not place the INQUA device immediately beside or on top of other devices.
- ⇒ Keep at a minimum distance of 30 cm from portable wireless communication devices (including accessories therefore, such as antenna cables or external antennas).

#### **NOTICE:Risk of damage to the compressor due to overheating and dust**

The compressor is fitted with a fan to ensure that the motor is cooled constantly. If the fan is not able to cool the motor sufficiently, the motor may overheat and be damaged.

- ⇒ Only operate the compressor in locations with low dust contamination and an adequate air supply. The ventilation slits of the compressor must always remain unobstructed.
- ⇒ Never cover the compressor while it is in use.
- ⇒ Do not run the compressor on carpets or in very dusty environments, e.g. under the bed, in workshops etc.
- ⇒ Observe the ambient conditions specified in section "Environmental conditions for operation", page 33.

The compressor is only intended for indoor use, it is not intended for outdoor use.

## **Information about electromagnetic compatibility:**

Medical electrical equipment is subject to special precautionary measures with regard to electromagnetic compatibility (EMC). Such equipment must only be installed and operated in accordance with EMC instructions.

Portable and mobile HF (high-frequency) communication devices can disrupt electrical medical equipment. Using accessories, converters and power cords other than those specified (with the exception of converters and power cords that the manufacturer of the medical electrical device sells as spare parts for internal components) can result in higher HF emission levels or lower the resistance to interference of the device.

The INQUA® Inhalator uses HF energy exclusively for its internal functions.

The HF emission from the compressor is very low. It is unlikely that it will interfere with other electronic equipment in the vicinity.

The electrical field generated during compressor operation can interfere with screens and monitors (flickering). For this reason, the compressor must be placed far enough away from screens/monitors so that image quality is not affected.

The INQUA® Inhalator is suitable for use in all buildings, including those used for accommodation and those that are connected directly to a public power supply which also supplies residential buildings.

Floors should be made from wood or concrete or covered with ceramic tiles. If the floor is covered with a synthetic material, the relative humidity must be at least 30%.

The quality of the supply voltage should be equal to that of a typical residential, business or hospital environment and must match the information on the type identification plate.

Magnetic fields in mains frequency should correspond to the values typically found in a residential, commercial or hospital environment. Technical data on electromagnetic compatibility (EMC information) is available in table format upon request from INQUA GmbH or on the internet at: [www.inqua.de/service/downloads/](http://www.inqua.de/service/downloads/)

**EMV\_Tabellen\_INQUA\_Inhalator.pdf.**

#### 4.3 Preparing for inhalation

##### **WARNING:Danger of electric shock**


Do not take hold of the compressor's power cord or power plug with wet or damp hands, and do not operate the compressor with wet or damp hands.

##### **CAUTION:Quality of treatment impaired by electromagnetic interference**

The use of third-party products can result in higher levels of electromagnetic emissions or lower the resistance of the INQUA device to electromagnetic interference.

- ⇒ Use only original spare parts and original accessories from INQUA.

Figure 3

- ⇒ Wind the power cord (4) up on the underside of the compressor.
  - ⇒ Place the compressor (1) on a flat, firm surface.
  - ⇒ Plug the power plug into the mains socket.
  - ⇒ Assemble the nebuliser (2) as described in the instructions for use of the nebuliser.
  - ⇒ Connect the nebuliser and the compressor to the connection tube (3) (as illustrated in Fig. 3). Ensure that the ends of the connection tube are connected firmly to the air connection (7) and the nebuliser.
  - ⇒ To fill the nebuliser, place it in the nebuliser holder (9) on the compressor.
- The INQUA® Inhalator is ready for use.

##### *Info:*

*The compressor  
can only be used in  
conjunction with the  
nebuliser recom-  
mended by INQUA.*

#### **4.4 Performing the inhalation**

##### **WARNING:Danger of death by electrocution from a faulty device!**

- If the device has a fault, a person may come into contact with live parts, resulting in the danger of electric shock, for example.
- ⇒ Switch the compressor off and disconnect the power adapter from the mains socket immediately if a fault is suspected (e.g., following a fall or if there is a smell of burning plastic).

##### **WARNING:Danger to health from inhaling unsuitable medications or incorrect dosages or incorrect use**

- ⇒ Follow the instructions in the package insert for the medication.
- ⇒ Before inhalation, read and comply with the instructions for use of the nebuliser.
- ⇒ Children and anyone who requires assistance must be supervised constantly by an adult during inhalation therapy. Otherwise, safe and effective treatment cannot be guaranteed.

##### **WARNING:Danger to health if connection tubings are mixed up**

Check carefully to ensure that the other end of the connection tubing connected to the compressor is attached to the nebuliser.

If multiple tube systems are used at the same time for different applications (e.g., infusion and inhalation), there is a risk that the wrong tube systems may be used by mistake, which in extreme cases may result in death.

##### **CAUTION:Danger of burning due to direct contact with the air connection during inhalation**

The ambient air, which is forwarded from the compressor to the connection tube heats the air connection on the compressor. As a result, the air outlet can become very hot. Therefore, avoid touching the air connection during an inhalation session.

##### **CAUTION:Quality of treatment impaired by electromagnetic interference**

If electrical devices are placed too close to each other, the function of both can be impaired due to electromagnetic interference.

- ⇒ If the INQUA device has to be placed immediately beside or on top of other devices for operation, all devices must be monitored to ensure that they are working properly.

*Info:*

*To ensure that the compressor functions perfectly, it must be at room temperature before it is put into operation.*

- ⇒ Fill the nebuliser as described in the instructions for use of the nebuliser.
- ⇒ Pull the nebuliser up and out of its holder.
- ⇒ Switch the compressor on at the ON/OFF switch (8). The nebuliser will now generate a clearly visible aerosol mist.
- ⇒ Perform the inhalation as described in the instructions for use of the accompanying nebuliser.

The inhalation session is complete when no more aerosol is produced by the nebuliser.


#### **4.5 Completing the inhalation**

*Info:*

*To ensure that it is completely disconnected from the mains power supply, the power plug must be unplugged from the mains socket.*

- ⇒ Place the nebuliser back in its holder.
- ⇒ Switch the compressor off.
- ⇒ Disconnect the power plug from the mains socket.
- ⇒ Clean the compressor (see Chapter "Cleaning, maintenance, repair", Section "Cleaning tasks", page 43).
- ⇒ Prepare the nebuliser and the mask for hygienic re-use as described in the instructions for use of the nebuliser.
- ⇒ Wind up the power cord in a clockwise direction on the bottom of the compressor.
- ⇒ When they have dried, stow the parts in the storage compartment of the compressor.

Figure 4

- ⇒ Keep the INQUA® Inhalator in a dry, dust-free place out of the reach of children (see Chapter "Miscellaneous", Section "Storage", page 46).

## 5 CLEANING, MAINTENANCE, REPAIR

Cleaning, maintenance and repairs to the device must not be carried out unless it is switched off and the power plug has been disconnected from the main socket.

### 5.1 Cleaning tasks

#### **WARNING: Danger of electric shock or malfunction caused by liquids getting into the compressor**

Liquids can conduct electric current. If liquid has got into the compressor, a current may then be present on the housing, causing a short circuit in the compressor.

- ⇒ Never spray any liquids onto the compressor or into the compressor's ventilation slits.
- ⇒ Wipe the compressor down with a damp cloth, not a wet one.
- ⇒ Do not immerse the compressor in liquids such as water or disinfectant solutions.
- ⇒ Always switch the compressor off and disconnect the power plug from the mains socket before working on the inhalation device.  
Always pull on the plug to disconnect it, not the cord.

If liquid has got into the compressor, do not use the compressor again until you have contacted the manufacturer (see back of instructions for use).

#### **CAUTION: Danger of infection from contaminated nebuliser/mask**

- ⇒ Read the chapter entitled "Hygienic re-use" in the instructions for use of the nebuliser you are using, and carry out the procedures for hygienic re-use described there.
- ⇒ Wipe the compressor housing and the storage compartment down with a dry or slightly damp cloth.
- ⇒ Carry out the steps for hygienic re-use for all other components as described in the instructions for use of the nebuliser.

**5.2 Maintenance** With the exception of the filter, the compressor contains no parts that are serviceable by the user.

If used twice daily for inhalation for a total of twice 20 minutes a day, the expected life cycle of the compressor is four years. If the compressor is still in use after four years, contact INQUA to arrange to have it inspected.

### Filter change

#### **CAUTION: Damage to the compressor or impaired treatment caused by unsuitable air filter**


If an air filter from a third party manufacturer is used, the compressor may be damaged or the inhalation treatment may be impaired.

⇒ Use only air filters that are recommended by INQUA for use with the INQUA® Inhalator.

The air filter in the compressor must be replaced after every 200 operating hours, but at least once a year or if it becomes visibly dirty, discoloured or damp.

To do this, proceed as follows:

Figure 5

#### **6 Air filter**

#### **11 Removal tool for changing the filter**

- ⇒ Fit the removal tool (11) firmly on the filter (6) in the compressor.
- ⇒ Turn the removal tool until you can pull the old filter out of the compressor.
- ⇒ Pull the old filter out of the removal tool and dispose of it.
- ⇒ Push the new filter as far as it will go into the filter socket with the opening facing forwards.

### Replacing the nebuliser

Note the information regarding material resistance and replaceability in the instructions for use of the nebuliser.

---

### **5.3 Repairs**

*Info:  
If the compressor  
is opened without  
authorisation, all guar-  
antees and claims  
under the warranty will  
be null and void.*

Repairs must not be carried out on the INQUA® Inhalator except by the manufacturer. If the device needs to be repaired, contact the manufacturer (see back of instructions for use).

---

### **5.4 Storage**

Keep the compressor and all accessories out of continuous direct sunlight.

Do not leave the compressor connected to the mains when storing because all electrical devices that are connected to mains electricity represent a potential hazard.

## **6 TROUBLESHOOTING**

Fault	Possible cause	Procedure
The compressor does not start up.	The power plug is not plugged into the mains socket properly.	Check that the power plug is seated correctly in the mains socket.
No aerosol is coming out of the nebuliser.	The nebuliser nozzle is blocked. The connection tubing is not connected properly.	Clean the nebuliser. Check that both ends of the tubing are seated properly.
	The connection tubing is leaking.	Replace the connection tubing.

Contact the manufacturer

- if the motor is unusually loud,
- if it switches off automatically after a few seconds,
- if no compressed air is generated,
- if the recommended procedure is unsuccessful.

## 7 MISCELLANEOUS

- 7.1 Transportation** For transporting inside a building, stow the parts of the INQUA® Inhalator in the compressor's storage compartment and hold the INQUA® Inhalator by the carrying handle.  
For transporting over longer distances, pack the INQUA® Inhalator in its original packaging. This is the best way to protect the product from damage.

<b>7.2 Storage</b>	If the INQUA® Inhalator is to be stored for a prolonged break in treatment, the following environmental conditions must be observed:
Minimum ambient temperature:	-25 °C (without monitoring of relative humidity)
Maximum ambient temperature:	+70 °C (with relative humidity of up to 93%, non-condensing)
Humidity:	0% – 93%
Atmospheric pressure:	500 hPa to 1,060 hPa

Store the compressor and accessories safely out of continuous direct sunlight.

Additional instructions for storing the nebuliser are included in the instructions for use of the nebuliser.

- 7.3 Terms and conditions of warranty** We guarantee that your product, if used according to the instructions, will be free from defects in material and workmanship caused by the manufacturing process for the indicated warranty period of 24 months (2 years) beginning on the date of initial purchase. The warranty period will expire after 12 months. This warranty applies in addition to the warranty obligation of your vendor. Your statutory rights with respect to your dealer in the event of defects are not limited by the warranty or any claim under the warranty. Your receipt serves as your proof of warranty and ownership.

### ***What does the warranty cover?***

If, exceptionally, a defect is discovered, we will at our discretion repair or replace the product, or refund the purchase price of the product. If it is replaced, the replacement product may either be the same model or a model that is at least comparably equipped.

Replacement or repair of the product shall not serve as the basis for a new warranty. All used devices or parts that are replaced become our property. Further claims are excluded. This applies particularly for any claims for compensatory damages. This disclaimer of warranty shall be ineffective in the event of injury to life, limb or health, in cases of wilful malfeasance and gross negligence, product liability and if substantive obligations under the warranty agreement are violated.

### ***The warranty shall be cancelled if***

- the product has been operated or used improperly with respect to the descriptions in the instructions for use
- damage is present that is attributable to the effects of water, fire, lightning, or the like
- the damage was caused by transporting the device incorrectly or a falling impact
- the product has been misused or not cared for correctly
- the serial number on the compressor has been changed, removed, or otherwise rendered illegible
- repairs, adaptations or modifications have been made to the product by unauthorised persons

Moreover, the warranty does not cover wearing parts, that is to say device parts that are exposed to normal wear. In the event of a claim under the warranty, please contact the manufacturer directly (see back of instructions for use).

The "warranty period" commences on the date of purchase. Please keep your receipt as your proof of purchase.

---

<b>7.4 Disposal</b>	This product falls within the scope of the European Council Directive on Waste Electrical and Electronic Equipment (WEEE <sup>1</sup> ) and is included in Product Category 8: Medical Devices for the purposes thereof. Accordingly, this product must not be disposed of with domestic waste. The disposal regulations prevailing in the respective member countries must be observed (e.g., disposal by local authorities or dealers). Materials recycling helps to reduce the consumption of raw materials and protect the environment.
<b>7.5 Information about RoHS</b>	This product satisfies the requirements of Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (RoHS). The product is therefore marked with the characters "CE0123".

## 8 SPARE PARTS

For spare parts, please contact the manufacturer (see back of instructions for use).

Information as of: July 2017

1. *Directive 2002/96/EC of the EUROPEAN PARLIAMENT AND THE EUROPEAN COUNCIL of January 27, 2003 on waste electrical and electronic equipment.*

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